



IHE Integration Statement			
Vendor	Product Name	Version	Date
Ivenix	Ivenix	1.0	9/24/14

This product implements all transactions required in the IHE Technical Framework to support the IHE Integration Profiles, Actors and Options listed below :

Integration Profile	Actor	Integration Profile Option
Alert Communication Management	Alarm Reporter	None
Consistent Time	Time Client	None
Device Enterprise Communication	Device Observation Reporter	None
Infusion Pump Event Communication	Device Observation Reporter	None

Internet address for vendors IHE information : <http://www.ivenix.com>

Links to Standards Conformance Statements for the Implementation

HL7	
DICOM	

Link to general information on IHE : www.ihe.net