

For Hospital IT

# Ivenix Infusion System

## IHE Integration Statement

Vendor: Fresenius Kabi  
Product Name: Ivenix Infusion System  
Version: 1.0  
Date: January 2018

**This product (the Ivenix Infusion System) implements all transactions required in the IHE technical framework to support the IHE integration profiles, actors, and options listed below:**

Integrated Profiles Implemented	Actors Implemented	Type
DEC - Device Enterprise Communication	Device Observation Reporter	Thorough
DEC - With ADT Demographics	Device Observation Reporter	Thorough
PIV - Point of Care Infusion Verification	Infusion Order Consumer	Thorough
ACM - Alert Communication Management	Alert Reporter	Thorough
IPEC - Infusion Pump Event Communication	Device Observation Reporter	Thorough
MEMLS - Medical Equipment Management Location Services	Location Observation Consumer	Thorough
Consistent Time	Time Consumer	Thorough

The internet address for Fresenius Kabi's IHE information is: <https://fresenius-kabi.us/Ivenix/IHE>

The links and standards conformance statements for HL7 and DICOM are not applicable.

Publication of this statement: July 2022

Link to general information on IHE for North America: <https://www.ihe.net>

Link to general information on IHE for Europe: <https://www.ihe-europe.net>

**IVENIX™**  
INFUSION SYSTEM