



## DATA SHEET

# Large Volume Pump (LVP)

Ivenix Infusion System

## Compatible Cleaning Solutions and Disinfectants

Proper care and maintenance of the Ivenix Large Volume Pump (LVP) will help ensure optimal performance and longevity of the device. Fresenius Kabi recommends that the LVP is cleaned and sanitized between patient use and prior to maintenance. The following is a list of recommended and compatible cleaning solutions and disinfectants.

### Cleaning Solutions

- EcoLab® Asepti-Zyme™ Multi-Enzymatic Detergent
- ENZOL® Enzymatic Detergent
- Metrex™ Metrizyme™ Dual-Enzymatic Detergent
- STERIS® Valsure® Alkaline Detergent
- STERIS Valsure Neutral Detergent

### Disinfectants

- Metrex CaviCide™ \*
- Clorox® Regular Bleach2 \*
- Clorox Healthcare® Hydrogen Peroxide Cleaner Disinfectant Wipe \*
- PDI® Sani-Cloth® Bleach Germicidal Disposable Wipe \*
- PDI Super Sani-Cloth Germicidal Disposable Wipe \*
- Wexford Innovations Ready To Use Thymo-Cide

### Note:

- Follow manufacturer's instructions on dilutions (if concentrate) and on recommended dwell times.
- Consult the Instructions For Use (IFU) for instruction and cautions for cleaning.
- Cleaning solutions and disinfectants effective against SARS-CoV-2 are indicated by an asterisk (\*).

For questions, contact [info@ivenix.com](mailto:info@ivenix.com).

Distributed in the U.S. by  
Fresenius Kabi USA, LLC  
Three Corporate Drive  
Lake Zurich, IL 60047  
Phone: 1.888.386.1300  
[www.fresenius-kabi.us/ivenix](http://www.fresenius-kabi.us/ivenix)

Manufactured for: Fresenius Kabi, AG  
61352 Bad Homburg, Germany

 **FRESENIUS  
KABI**  
caring for life